



LITHIUM-ION BATTERY

- Faster charging time
- Sustained power performance
- Lighter weight
- More usable energy
- Longer life
- Smart battery protection

Power Chairs	Battery Size	Total Watt
Jazzy® Passport* (Travel Battery)	25.2 V 10.4 Ah	262Wh
Jazzy® Passport**	25.2 V 18.2 Ah	458Wh
Jazzy® EVO 613 Li**	25.6 V 32.4 Ah	829Wh
Jazzy® Carbon*	24 V 12Ah	288Wh
Scooter/Travel	Battery Size	Total Watt
i-Go™	25.2 V 9.6 Ah	242 Wh
iRide® 2*	37 V 5.2 Ah	192Wh
Go Go® Endurance Li 8Ah*	25.6 V 8 Ah	204Wh
Go Go® Endurance Li 16Ah Maxx**	25.6 V 16 Ah	409Wh
Go Go® Folding Scooter*	24 V 11.5 Ah	276Wh

Note : Veuillez consulter le 49 CFR 175.10(a)(17) pour toute exigence ou condition supplémentaire susceptible de s'appliquer.

*Remarque : le transport de cet article à bord d'un avion commercial est soumis aux réglementations locales et à la politique de la compagnie aérienne. Veuillez consulter votre compagnie aérienne avant votre vol pour vous assurer que votre fauteuil ou scooter et ses batteries lithium-ion sont autorisés.

**Ces batteries ne sont pas conformes aux normes des compagnies aériennes.